

MRI

The ExoDx™ Prostate Test & mpMRI - A Complementary Approach

Contents

Introduction	2
Abstract	3
Clinical Challenge	3
Multiparametric Magnetic Resonance Imaging (mpMRI)	3
Areas for Enhancement	4
mpMRI Performance	6
ExoDx™ Prostate	8
mpMRI and ExoDx™ Prostate	10
Conclusions	12
References	14

Author: Dr. Jason M. Alter is the Head of Scientific and Clinical Affairs for ExosomeDx, a Bio-Techne brand. He has extensive experience in prognostic and predictive testing in Urologic Oncology and has played key roles in the development of prostate cancer assays (biopsy and post-radical prostatectomy) at multiple companies. Dr. Alter has published on genomic and non-genomic risk assessment methods for prostate cancer, and he is very familiar with available molecular and genomic testing for the disease. Dr. Alter has a dual B.A. in Biology and History from Alfred University, an M.S. degree in Immuno-Parasitology from Texas A&M University and a Ph.D. in Molecular Biology from Binghamton University. He did a postdoctoral fellowship at Schering Plough Pharmaceuticals.



Introduction

One of the most widely adopted risk assessment tools employed in prostate cancer (PCa) management is multiparametric magnetic resonance imaging (mpMRI). The EAU, NCCN, and AUA guidelines call for mpMRI utilization in multiple places in the prostate cancer management pathway. Protocols have been implemented, although some have expressed caution about moving too swiftly.^{1-3, 80-82} The NCCN Guidelines for Early Detection of Prostate Cancer currently recommend that mpMRI precede biopsy (naïve or prior negative biopsy).

Risk assessment tools including mpMRI all have strengths and limitations. mpMRI limitations include well-known dependencies on reader expertise and therefore variable interpretation, tumor size, multifocality, tissue architecture, process standardization, ethnic and racial availability disparity, potential water resource contamination, and, of course, cost.⁴⁻⁶ Patient management is best informed when complementary risk assessment methods are appropriately combined to mitigate the limitations of each approach (FIGURE 1).

In this whitepaper, we review the potential strengths and limitations of mpMRI and how biomarkers, specifically the ExoDx Prostate Intelliscore Test (also known as EPI), might be considered, in combination with imaging, to provide the most informed decision-making.^{7,8}

Abstract




Figure 1 Potential clinical strategies for combining the ExoDx Prostate Test and multiparametric magnetic resonance imaging (mpMRI). Using the ExoDx Test either before or after mpMRI are possible approaches for layered risk assessment and more informed decision-making.^{79, 82, 83}

Clinical Challenge

In 2022, 268,490 newly diagnosed prostate cancer cases are projected, and ~34,500 men are expected to die from the disease.⁹ However, since the United States Preventative Services Task Force's (USPSTF) decision to alter PSA recommendations, the rate of distant metastasis has increased,¹⁰⁻¹² and there is concern that not doing sufficient biopsies might miss high-grade prostate cancer (HGPCA) (\geq GG2). Although prostate biopsy should be employed to find HGPCA, guidelines and much of the medical community look to minimize the detection of low-risk PCa (GG1). Most of these patients will be fine for many years, and if diagnosed with GG1, most will be directed to active surveillance (AS).^{1,3}

The clinical question before us is how do we enable the best clinical balance, providing biopsies to men at higher risk for HGPCA while deferring biopsy in men with low risk for HGPCA?

Multiparametric Magnetic Resonance Imaging (mpMRI)

mpMRI is a powerful technology that provides additional insight into which tumors may be clinically significant and where they might reside. Many studies demonstrate the high sensitivity of mpMRI detection for HGPCA, ranging from 82% to 100%, while specificity ranges from 35% to 100%.¹³⁻²¹ Tumor detection (sensitivity) depends on the type of biopsy performed (TRUS, targeted fusion, targeted visual, etc.). Though still debated, it is widely accepted that mpMRI imaging with biopsy provides either non-inferior or superior detection over TRUS biopsy.^{14, 16, 19, 22-24, 25-35} NCCN guidelines highlight the value of mpMRI and the debate about the appropriate mix of image-guided with standard biopsy for tumor detection.³



Areas for Enhancement

There are acknowledged variations in mpMRI reporting, targeted biopsy methodology, and associated protocols.⁵¹

1. Reader Variability/Subjectivity

One of the most widely understood mpMRI limitations is reader-dependent variation. In some studies using experienced or 'expert' readers, there are claims of high concordance (78%).²³ Notably, Rosenkrantz *et al* used experienced readers and observed moderate interobserver variability for \geq PIRADS 4 with better agreement in the posterior zone than the transitional zone.⁵² Greer *et al* used whole-mount radical prostatectomy (RP) pathology to note that highly experienced readers did better than moderately skilled readers only for \geq PIRADS 4 lesions - not for \geq PIRADS 3 tumors. They suggest that less experienced readers were likelier to score the lesions as PIRADS 3 due to the ambiguity in PIRADS 3 and PIRADS 4 lesions.⁴ Standardization is critical to consistent, agreed-upon MRI performance metrics. In a large multinational study spanning 26 sites, the positive predictive value (PPV), the likelihood that a positive result for HGPCA (\geq GG2) is truly positive, varied widely from 15%, 39%, and 72% for PIRADS 3, PIRADS 4 and PIRADS 5 respectively.⁶ Variation in PPV was likely the result of reader variation, poor targeting, and inconsistent disease prevalence across sites.

2. Multifocality/Disease Heterogeneity

Prostate cancer is often multifocal with different foci displaying both Gleason and genomic heterogeneity.^{5,46,47} Unfortunately, mpMRI does not do as good a job detecting multifocal tumors as solitary lesions; multifocality increases the probability of mpMRI-missed tumors^{5,34,39} (FIGURE 3). A literature review reveals that 11%-13% of men had HGPCA at a different prostate location than that captured by mpMRI; this raises particular concerns for focal therapies.⁴⁸ Stabile *et al* followed fusion biopsy with TRUS biopsy and found 30% HGPCA outside of the index lesion: the Gleason score was equal to or higher than that found in the index lesions. They observed that the likelihood of finding HGPCA outside the MRI-located index lesion grew as the PIRADS score increased: 10% probability for PIRADS 2 and 70% probability for PIRADS 5.³⁴ Studies that incorporate whole-mount radical prostatectomy (RP) are the ultimate pathological ground truth and provide insight into the impact of tumor size and multifocality on mpMRI detection.^{5,39} Bonekamp *et al* observed that even though MRI-detected HGPCA lesions were detected in 97% of men, there were additional \geq PIRADS 3 lesions that were missed in 60% of these cases.²⁵

Tumors Detected on mpMRI (%)




Figure 2 The impact of tumor size on mpMRI tumor detection. Bubble size reflects the percentage (%) of tumors detected by mpMRI: larger bubbles indicate a higher detection rate than smaller bubbles.

- *Sizes based on solitary tumors
- ⁸Sizes based on multifocal tumors
- Green bubbles are mixed Gleason grade groups while blue bubbles are exclusively high-grade disease (≥GG2)^{39-41, 44, 49}

3. Tumor Size

Tumor size is a critical factor for mpMRI detection. Studies employing whole-mount radical prostatectomy demonstrate the relationship between tumor size, multifocality, and tumor detection.^{5,39,40,41} Although tumors greater than 1 cm can be missed by mpMRI, tumor detection decreases as tumor size decreases.^{39-42,44} Also, small tumors are understood to harbor high-grade disease.^{5,42,44}

4. Invisible Tumors

Up to 55% of all prostate cancer and 35% of HGPCA tumors are not visible on mpMRI.^{5,28,43} In addition to size, grade, and multifocality, specific tissue architecture (cribriform) reduces tumor visibility.⁴⁹ Moreover, gene expression appears to be involved with poor cellular organization in tumor visibility, at least to some degree. Interestingly, visible and invisible tumors have similar outcomes: biochemical recurrence, distant metastasis, and prostate cancer-specific mortality and do not impact genomic test scores.^{50,47}

5. PIRADS Variation

PIRAD scores are risk assessment categories that share the limitations of all risk assessment tools regardless of type - biomarker, nomogram, or tissue-based test.

A PIRADS category provides a probability - not a guarantee - of the result to be found at biopsy. Studies that have biopsied all PIRADS results have displayed the diverse range of HGPCA in each category, low-risk cancer (GG1), and benign tissue that makes up varying degrees of each PIRADS category.^{3, 34, 57, 74, 76}

6. False Positives

Many factors impact appropriate MRI-specific tumor detection. Conditions such as hyperplasia, inflammation, fibrosis, prostatitis, and high-grade prostatic intraepithelial neoplasia (HGPIN) can cause false positive MRI readings.³⁶

7. False Negatives

Tumor size, multifocality and tissue architecture can all contribute to false negatives. But, tumor location is also important to mpMRI detection. Studies find lesions in all zones - peripheral, transition, and anterior regions.^{13,25,28,37,38} Serrao *et al*, demonstrated that 50% of mpMRI false negatives were located in the anterior region. In contrast, Kido *et al*, utilized radical prostatectomy to note that 53% and 47% of MRI undetected tumors were in the peripheral and transition zones, respectively.^{37,38}

mpMRI Detection of Multifocal Prostate Cancer




Figure 3 Studies investigated mpMRI ability to detect multifocal disease. Although mpMRI misses the index lesion 20%-30% of the time, non-index lesions are missed much more frequently (63%-79%).

mpMRI Performance

Comparing mpMRI performance across studies is difficult because of study variation, which includes everything from different equipment (e.g., magnet strength) to consistent incorporation of the three (3) individual phases (T2 weighted, diffusion, and DCE), as well as reader experience. Study design is as guilty, if not more so, than equipment standardization for much of the difficulty in assessing and comparing MRI metrics. Biopsy methodology varies greatly, as does which samples are evaluated; studies often mix biopsy naive cases with prior negative biopsies and/or samples from men on active surveillance.^{25, 53-55} These 'mixed use' inclusion criteria impact disease prevalence which affects positive and negative predictive values.⁵³

Even the definition of HGPCA is not consistent across mpMRI studies. While many use ISUP grade group 2 or higher (\geq GG2), other definitions also include core length or percentage of cores positive, or ISUP grade \geq GG3.^{20, 56-61} Mortezaei *et al*, illustrates how HGPCA definitions and biopsy methodology impact HGPCA incidence.⁶²

mpMRI Negative Predictive Value (NPV)

As mentioned, mpMRI seems better at finding larger, higher-grade solitary tumors than multifocal tumors.^{5,39-41,60} Positive predictive values (PPVs), or the likelihood that the positive MRI will find HGPCA, vary but trend upward as PIRAD scores increase.^{23,51,58}

When it comes to negative or low-risk MRI results (PIRADS $<$ 3), the range of negative predictive values (NPV) is wide.^{20,25,43,51,53-62} Some data suggests 'good' mpMRI

NPV on a per patient basis but not on a per lesion basis.⁵⁷ Arguably, NPV is more concerning than PPV because negative mpMRI results are used to avoid biopsy,²⁸ and concerning HGPCA (\geq GG2) is often found after a negative or low-risk mpMRI.^{21,51,56} Chung *et al*, examined 'invisible' tumors (PIRADS $<$ 3) with both biopsy and radical prostatectomy (RP). They found that, at biopsy, 24% and 6.6% were \geq GG2 and \geq GG4 respectively.⁶³ One study even found that 17% of men with 'normal' mpMRI readings had palpable disease.⁵¹ Kinnaird *et al* examined men initially classified as mpMRI negative at 2-4 years post-imaging and found that the false negative rate was 23% for HGPCA.⁶⁴ Guidelines reflect published data emphasizing that, by itself, a negative mpMRI does not omit the possibility of cancer and clinicians should consider biomarkers when looking to defer a biopsy in a patient with a negative mpMRI.³

Early studies presented NPV claims based on small numbers of cases, muddled mixed biopsy sample types (biopsy naive, prior negative, active surveillance) without sub-analysis, varied disease prevalence, or a combination that impacted NPV.^{15,65} Moldovan *et al* did a thorough review of studies documenting how clinical heterogeneity and disease prevalence impacted NPV results for HGPCA that ranged from 64% to 88%. Therefore, we only included studies that were either published after the Moldovan *et al* review or were not cited in their study.⁶⁵

Studies continue to mix and match different biopsy sample types resulting in different disease prevalence that impacts PPV and NPV.^{25,53-55,61} Several also use alternate definitions of HGPCA that include core number or percent of core positive.^{20,56,60} Fortunately, some studies conduct sub-analysis to examine metrics separately based on biopsy type^{54,56,62,66} (FIGURE 4).

Studies that utilize a more holistic assessment of the prostate, employing template biopsy mapping, saturation biopsy, or whole mount radical prostatectomy, offer the most comprehensive pathologic ground truth for assessing performance metrics.^{5,56,57,60,67} One of the most comprehensive studies to provide a detailed assessment of mpMRI metrics was the PROMIS study.⁵⁶ The PROMIS study did many things well; high-quality and standardized MRI, in depth reporting, dedicated and experienced urologic radiologists, centralized reader training, as well as high-quality targeted biopsies. PROMIS researchers performed template mapping biopsy every 5 mm and a TRUS biopsy on every man: an approach that provides a much better understanding of the prostate than TRUS or fusion biopsy. One potential criticism of the PROMIS study was the use of 1.5 Tesla magnets instead of the 3 Tesla magnets most currently use. However, much of the essential mpMRI

literature has utilized a mix of magnet strengths, including PRECISION, MRI-FIRST, STHLM3-MRI (NCT03377881).^{22,23,33} Furthermore, some researchers have found no difference in PPV or NPV when comparing 1.5T or 3T generated data.^{6,20}

The PROMIS study had three definitions of clinically significant cancer: (1) \geq GG3 or cancer core length \geq 6 mm, (2) \geq GG2 or cancer core length \geq 4 mm, and (3) \geq GG2. Each include different performance metrics. The primary definition (\geq GG3 or cancer core length \geq 6 mm) resulted in a PPV of 51% and NPV of 89%. But \geq GG2 is widely utilized as the definition of HGPCA in clinical practice, and PROMIS observed that for \geq GG2, the PPV was 65% and NPV was 76%.

Other studies have utilized mapping biopsies and demonstrated similar NPVs for \geq GG2 or related definitions. Simmons *et al* noted an NPV of 68.6% for \geq GG2 and/or tumor length of \geq 4 mm.⁶⁰ Mortezaei *et al* used template saturation biopsy to measure mpMRI performance metrics

for \geq GG2 and noted overall NPVs of 74.2% and 68.5% for saturation biopsy and targeted fusion biopsy, respectively. They also did subgroup analysis by prostate sample type (naïve or prior negative biopsy) resulting in NPVs that were changed due to the impact of disease prevalence on NPV and PPV.⁶² In a retrospective analysis, Hogan recently observed an NPV of 95% for HGPCA and took a higher number of cores (median=26), but they did not utilize saturation biopsy. They hypothesize their improved measurements are because they employed one highly experienced, centralized reader.⁵⁹

mpMRI Cost

mpMRI is used widely in the United Kingdom, but mpMRI usage varies significantly from country to country in Europe; surveys show that cost is the main reason mpMRI is often not employed.⁶⁸ mpMRI is an expensive procedure with a median cost of \$4396 (IQR \$2,784-\$7,127) for MRI-guided biopsy with imaging contributing significantly (median of \$1,704) to that total. If anesthesia is utilized, the total cost increases to \$5,832.⁶⁹

mpMRI NPV




Figure 4 References reviewed in Moldavan *et al*.⁶⁴ are excluded except for the PROMIS study.⁵⁶ All NPVs reported for HGPCA are defined as \geq GG2 unless a specific notation is included. Ball location on the Y axis indicates NPV while ball size indicates cohort size.

*NPV for \geq GG2 with \geq 10% of cores positive or \geq 5mm in any core.

**One (1) skilled reader

***55% of cohort with initial biopsy, 30% with prior negative biopsy and 16% on active surveillance.

†81% of cohort with initial biopsy, 19% with prior negative biopsy.

‡39% of cohort with initial biopsy, 61% with prior negative biopsy.

ExoDx™ Prostate Test

Exosomes are very small (30-150 nanometers (nm)) vesicles containing DNA, RNA, and proteins and are produced and secreted in large quantities from every cell in the body. In healthy individuals, billions of exosomes exist per milliliter (mL) of plasma.^{70,71} In addition to blood, exosomes are found in all bodily fluids (urine, cerebrospinal, breast milk, saliva, etc.). Exosome signaling has a critical advantage over classical, soluble signaling pathways as exosome signaling does not decrease over distance. Integral to the development and propagation of multiple diseases, exosomes are involved in normal physiology and cancer alike.^{72,73}

The ExoDx Prostate Test is a non-invasive urine assay that does not require a digital rectal exam (DRE) and analyzes three exosome-located genomic markers that are associated with high-grade prostate cancer (HGPCA) defined as \geq GG2. The ExoDx Prostate Test is validated for use in men with initial or repeat biopsy (NCT02702856 and NCT03031418) who are 50 years or older and have PSAs between 2-10ng/mL.^{74,75} ExoDx Prostate is a continuous risk assessment tool that results in a risk score: the risk of high-grade prostate cancer (\geq GG2) increases as the score increases.

ExoDx Prostate Test helps make more informed decisions on which patients may likely benefit – or not – from a prostate biopsy. In a randomized controlled two-arm trial (level 1 evidence) – one arm used standard of care (control arm) to guide patient management. In contrast, in the other arm, clinicians received and integrated the ExoDx Prostate Test results into their shared decision-making with patients. This study demonstrated that the ExoDx Prostate Test was able to defer biopsy in men unlikely to find HGPCA, suggest which men were at higher risk and should have a biopsy, and consequently detect 30% more HGPCA in the ExoDx Prostate Test arm than in the standard of care arm.⁷⁶


The ExoDx Prostate Test Performance

Genomic assessment using the ExoDx Prostate Test is significant compared to any individual clinical feature or combination of features. In validation studies, the ExoDx Prostate Test was compared to an optimized standard of care model (not found in clinical practice), clinical risk calculators such as PCPT and ERSPC, as well as clinical features (PSA), and was significantly more accurate than either optimized standard of care clinical models, risk calculators or PSA. Adding the standard of care information to the ExoDx Prostate genomics did not improve test performance – the genomic information was significantly more accurate as measured by the area under the curve (AUC).


The ExoDx Prostate Test has validated high NPVs of 91.3% for \geq GG2 and 97% for \geq GG3, while the PROMIS study demonstrated 76% NPV for \geq GG2. It's important to note that the ExoDx Prostate Test, like most tests, was validated on biopsy which is the 'gold standard.' The PROMIS NPV was based on a template mapping biopsy which is more complete. However, the ExoDx Prostate Test was assessed in some patients that subsequently had RP, and we found that the ExoDx Prostate Test PPV improved while its NPV stayed the same.⁷⁶ These results must be investigated on a larger scale. Until a direct head-to-head comparison of mpMRI and the ExoDx Prostate Test in the same cohort is conducted, we will not know the comparative performance metrics. Until that study is performed, we will look at the best data available: the ExoDx Prostate Test validation studies and the PROMIS trial.

Logistics and Convenience

In comparison to mpMRI, the ExoDx Prostate Test is straightforward to use. The ExoDx Prostate Test only requires a urine sample and does not require a digital rectal exam (DRE). In fact, due to the ExoDx Prostate Test's ease of use, a home sample collection kit was developed and deployed early during the COVID-19 pandemic. The home collection kit will continue to be a boon to clinicians and patients as prostate biopsy is forecast to be the outpatient procedure most often to burden urologists.



ExoDx™ Prostate Test



MRI and ExoDx™ Prostate

Combining risk assessment methods - such as mpMRI and biomarkers - (when appropriately applied) can better inform clinical decision-making. Published data demonstrates that combining different types of biomarkers, be it PSA density, risk calculators, or genomic testing, can provide enhanced assessment when appropriately layered with mpMRI in a clinical pathway.^{53, 58, 61, 67}

Until mpMRI and biomarkers are incorporated into well-designed RCTs, we will not truly understand the most beneficial way to integrate the two technologies. Nevertheless, based on published information, we can hypothesize how biomarkers, specifically the ExoDx Prostate Test, may be combined with mpMRI to answer the question. de la Calle *et al* performed an observational analysis of men who had both mpMRI and several different biomarkers including the ExoDx Prostate Test. This work, independent of Exosome Diagnostics involvement or sponsorship, was originally published in an abstract form and eventually as a published article.^{7,79}

The authors developed a series of algorithms that provide perspective on the possible results of combining the ExoDx Prostate Test and mpMRI in various clinical paths.⁷ Recreated (FIGURE 5) de la Calle *et al* highlights that the ExoDx Prostate Test use alone might avoid a significant percentage of biopsies while only deferring <5% of HGPCA detection. Moreover, employing the ExoDx Prostate Test and mpMRI together (algorithm #4) avoids both likely unneeded mpMRIs and biopsies without delaying the detection of any additional HGPCA.⁷

Preliminary work at Exosome Diagnostics investigating the ExoDx Prostate Test and mpMRI showed an association between rising ExoDx Prostate Test scores and PIRAD scores, (FIGURE 6). Our research also demonstrates the potential benefits of modeling the ExoDx Prostate Test and mpMRI together (FIGURE 7). Based on a retrospective analysis of 93 men aged ≥50 years with PSA 2-10 ng/ml under consideration for a prostate biopsy that received an MRI, an ExoDx Prostate Test, and a subsequent biopsy we determined that the ExoDx Prostate Test and MRI provide independent sources of information and in combination perform better than individually.

Possible Integration of ExoDx Prostate Test, mpMRI, and PSAD




Figure 5 Figure recreated from de la Calle *et al* 2020. In this study, independent of Exosome Diagnostics involvement or sponsorship, synergy between the ExoDx Prostate Test and mpMRI was shown by de la Calle *et al*. The study compared PSAD alone, ExoDx alone, mpMRI alone, PSAD + mpMRI, ExoDx + mpMRI + PSAD and an ExoDx + MRI algorithm (defer mpMRI and Bx below the ExoDx cutpoint (15.6), but have mpMRI if the ExoDx result is between 15.9 and 19 but Bx only if the MRI is positive. If ExoDx>19, obtain MRI and Bx regardless of MRI results). Elevated PSA was considered <20ng/mL.

ExoDx Prostate Test and mpMRI Association




Figure 6 Do ExoDx and MRI provide independent or the same (correlated) information? ExoDx and MRI are neither perfectly correlated nor orthogonal. Preliminary data indicates that the prevalence of both higher Gleason grade groups and higher ExoDx Scores increase with higher PIRADS Scores indicating a degree of correlation.

Modeling Potential ExoDx Prostate Test/mpMRI Combination




Figure 7 Preliminary modeling of combining ExoDx and mpMRI suggests superior area under the curve (AUC), with AUCs of 0.87 and 0.80 for the PIRADS +ExoDx model or PIRADS alone respectively. By DeLong's test, PIRADS+ExoDx performs better than MRI ($p=0.023$). This data supports the complementary value of a combined biomarker/mpMRI model.


Conclusions

All risk assessment methods, including mpMRI or biomarkers, have strengths and limitations. The key to more informed clinical decision-making is to really understand each technology's limitations and consider the appropriate integration of complementary risk assessment methods. There is a vital need to understand better how mpMRI and various biomarkers can provide more value to clinical practice; indeed, the newly formed ReIMAGINE Consortium was explicitly born to develop risk assessment tools that can examine the benefits of combining mpMRI with biomarkers.⁸

Biomarkers, in particular the ExoDx Prostate Test, have a complementary role with mpMRI. Combining the ExoDx Prostate Test with mpMRI has potential benefits for maximizing detection of HGPCA while minimizing HGPCA that may be missed by either method alone.

Author

Jason M. Alter, PhD
Head of Scientific & Clinical Affairs
Exosome Diagnostics, a Bio-Techne brand
Email: jason.alter@bio-techne.com



References

1. American Urological Association 2015. Early Detection of Prostate Cancer. AUA Guideline J Urol 2013 [cited 2021];190:419.
2. Mottet N, van den Bergh RC, Biers E et al. EAU Guideline Prostate. Cancer, Limited Update March 2022. Available at: <http://uroweb.org/guideline/prostate-cancer/>. Accessed July 2022.
3. Prostate Cancer Early Detection. National Comprehensive Cancer Network (NCCN) version 1.2023 – January 9, 2023.
4. Greer MD, Brown AM, Shih JH, et al. Accuracy and agreement of PIRADSv2 for prostate cancer mpMRI: a multireader study. J Magn Reson Imaging 2017;45(2):579-585.
5. Johnson DC, Raman SS, Mirak SA, et al. Detection of individual prostate cancer foci via multiparametric magnetic resonance imaging. Eur Urol 2019;75(5):712-720.
6. Westphalen AC, McCulloch CE, Anaokar JM, et al. Variability of the positive predictive value of PI-RADS for prostate MRI across 26 centers: experience of the Society of Abdominal Radiology Prostate Cancer Disease-focused panel. Radiology 2020;296(1):76-84.
7. de la Calle C, Fasulo V, Cowan JE, et al. Clinical Utility of 4Kscore, ExosomeDx™ and Magnetic Resonance Imaging for the Early Detection of High Grade Prostate Cancer. J Urol 2020;205:452-460.
8. Marsden T, McCartan N, Brown L, et al. The ReIMAGINE prostate cancer risk study protocol: A prospective cohort study in men with a suspicion of prostate cancer who are referred onto an MRI-based diagnostic pathway with donation of tissue, blood and urine for biomarker analyses. PLOS ONE 2022;2-18. doi.org/10.1371/journal.pone.0259672.
9. Siegel RL, Miller KD, Fuchs HE, and Jemal A. Cancer Statistics 2023. CA CANCER J CLIN 2023;73(1):17-48.
10. Jemal A, Culp MB and Fedewa SA. Prostate Cancer Incidence 5 Years After US Preventive Services Task Force Recommendations Against Screening. J Natl Cancer Inst 2021;113(1):64-71.
11. US Preventive Services Task Force. Screening for prostate cancer: US Preventive Services Task Force recommendation statement. JAMA 2018;319(18):1901-1913. doi:10.1001/jama.2018.3710.
12. Weiner AB, Matulewicz RS, Eggner SE and Schaeffer EM. Increasing incidence of metastatic prostate cancer in the United States (2004–2013). Prostate Cancer and Prostate Disease 2016; doi:10.1038/pcan.2016.30.
13. Abd-Alazeez M, Ahmed HU, Arya M et al. The accuracy of multiparametric MRI in men with negative biopsy and elevated PSA level – Can it rule out clinically significant prostate cancer? Urol Oncol Sem Orig Invest 2014;32:45e.17-22.
14. Freifeld Y, Xi Y, Passoni N, et al. Optimal sampling scheme in men with abnormal multiparametric MRI undergoing MRI-TRUS fusion prostate biopsy. Urol Oncol Sem and Orig Invest 2019;37:57-62.
15. Girometti R, Bazzocchi M, Como G, et al. Negative predictive value for cancer in patients with “gray-zone” PSA level and prior negative biopsy: preliminary results with multiparametric 3.0 Tesla MR. J Magn Reson Imaging 2012;36:943-950.
16. Haffner J, Lemaître L, Puech P, et al. Role of magnetic resonance imaging before initial biopsy: comparison of magnetic resonance imaging-targeted and systematic biopsy for significant prostate cancer detection. BJU Int. 2011;108(8 Pt 2):E171-E178.
17. Lamb BW, Shen Tan W, Rehman A, et al. Is Prebiopsy MRI Good Enough to Avoid Prostate Biopsy? A Cohort Study Over a 1-Year Period. Clin Genito Cancer 2015;13(6):512-517.
18. Lo G, Burton KR, Haider MA, et al. Negative predictive value of prostate multiparametric magnetic resonance imaging among men with negative prostate biopsy and elevated prostate specific antigen: a clinical outcome retrospective cohort study. J Urol 2019;202:1159-1165.
19. Siddiqui MM, Rais-Bahrami S, Turkbey B, et al. Comparison of MR/ultrasound fusion-guided biopsy with ultrasound-guided biopsy for the diagnosis of prostate cancer. JAMA 2015;313(4):390-397.
20. Thompson JE, van Leeuwen PJ, Moses D, et al. The Diagnostic Performance of Multiparametric Magnetic Resonance Imaging to Detect Significant Prostate Cancer. J Urol 2016;195:1428-1435.
21. Washino S, Okochi T, Saito K, et al. Combination of prostate imaging reporting and data system (PI-RADS) score and prostate-specific antigen (PSA) density predicts biopsy outcome in prostate biopsy naïve patients. BJU Int 2017; 119:225-33.
22. Eklund M, Jaderling F, Discacciati A, et al. MRI-targeted or standard biopsy in prostate cancer screening. NEJM 2021; 385:908-920.
23. Kasivisvanathan V, Rannikko AS, Borghi M, et al. MRI-targeted or standard biopsy for prostate-cancer diagnosis. NEJM 2018;378:1767-77.
24. Pagniez MA, Kasivisvanathan V, Puech P, et al. Predictive factors of missed clinically significant prostate cancers in men with negative magnetic resonance imaging: a systematic review and meta-analysis. J Urol 2020;204:24-32.
25. Bonekamp D, Schelb P, Wiesenfarth M, et al. Histopathological to multiparametric MRI spatial mapping of extended systematic and MR/TRUS-fusion-targeted biopsy of the prostate. Eur Radiol 2019;29:1820-1830.
26. Borkowetz A, Platzek I, Toma M, et al. Direct comparison of multiparametric magnetic resonance imaging (MRI) results with final histopathology in patients with proven prostate cancer in MRI/ultrasonography-fusion biopsy. BJU Int. 2016;118:213-220.
27. Dell'Oglio P, Stabile A, Soligo M, et al. There is no way to avoid systematic prostate biopsies in addition to multiparametric magnetic resonance imaging targeted biopsies. Eur Urol Oncol 2020;3: 112-8.
28. Filson CP, Natarajan S, Margolis DJA, et al. Prostate cancer detection with magnetic resonance ultrasound fusion biopsy: the role of systematic and targeted biopsies. Cancer 2016;122 (6):884-892.
29. Klotz L, Chin J, Black PC, et al. Comparison of multiparametric magnetic resonance imaging-targeted biopsy with systematic transrectal ultrasonography biopsy for biopsy-naïve men at risk for prostate cancer: a phase 3 randomized clinical trial. JAMA Oncol 2021; 7: 534-42.
30. Mannaerts CK, Kajtazovic A, Lodeizen OAP et al. The added value of systematic biopsy in men with suspicion of prostate cancer undergoing multiparametric MRI-targeted biopsy. Urol Oncol 2019;37: 298e1-e9.
31. Miah S, Hosking-Jervis F, Connor MJ et al. A multicentre analysis of the detection of clinically significant prostate cancer following transperineal image-fusion targeted and nontargeted systematic prostate biopsy in men at risk. Eur Urol 2020;3:262-269.
32. Ploussard G, Beauval JB, Lesourd M, et al. Added value of concomitant systematic and fusion targeted biopsies for grade group prediction based on radical prostatectomy final pathology on positive magnetic resonance imaging. J Urol 2019;202(6):1182-1187.
33. Rouviere O, Puech P, Renard-Penna R et al. Use of prostate systematic and targeted biopsy on the basis of multiparametric MRI in biopsy-naïve patients (MRI-FIRST): a prospective, multicentre, paired diagnostic study. Lancet Oncol 2019; 20:100.
34. Stabile A, Dell'Oglio P, De Cobelli F, et al. Association between prostate imaging reporting and data system (PI-RADS) score for the index lesion and multifocal, clinically significant prostate cancer. Eur Urol 2018;1:29-36.
35. Wang NN, Teslovich NC, Fan RE, et al. Applying the PRECISION approach in biopsy naïve and previously negative prostate biopsy patients. Urol Oncol Sem Invest 2019;37:530.e19-e24.
36. Gold S, Hale G, Bloom J, et al. Follow-up of negative MRI-targeted prostate biopsies: when are we missing cancer? World J Urol 2019;37:235-241.
37. Kido A, Tamada T, Kanomata N, et al. Multidimensional analysis of clinicopathological characteristics of false-negative clinically significant prostate cancers on multiparametric MRI of the prostate in Japanese men. Jap J Radiol 2019;37:154-164.
38. Serrao EM, Barrett T, Wadhwa K, et al. Investigating the ability of multiparametric MRI to exclude significant prostate cancer prior to transperineal biopsy. CUAJ 2015;9:11-12.
39. Le JD, Tan N, Shkolyar E, et al. Multifocality and prostate cancer detection by multiparametric magnetic resonance imaging: Correlation with whole-mount histopathology. Eur Urol 2015;67:569-576.
40. Wang M, Janaki N, Buzzy C, et al. Whole mount histopathological correlation with prostate MRI in Grade I and II prostatectomy patients. Int Urol and Nephrol 2019;51:425-434.
41. Wibulpolprasert P, Raman S, Margolis D, et al. Detection and Localization of Prostate Cancer at 3-T Multiparametric MRI Using PI-RADS Segmentation. AJR 2019;W122-W131.
42. Park MY, Park KJ, Lim B, et al. Comparison of biopsy strategies for prostate biopsy according to lesion size and PSA density in MRI-directed biopsy pathway. Ad Radiol. 2020;45:4166-4177.

43. Radtke JP, Kuru TH, Boxler S, *et al.* Comparative analysis of transperineal template saturation prostate biopsy versus magnetic resonance imaging targeted biopsy with magnetic resonance imaging ultrasound fusion guidance. *J Urol* 2015;193:87.
44. Arslan A, Karaarslan E, Guner AL, *et al.* Comparison of MRI, PSMA PET/CT, and fusion PSMA PET/MRI for detection of clinically significant prostate cancer. *J Comput Assist Tomogr* 2021;45(2):210-217.
45. Pooli A, Johnson DC, Shirk J, *et al.* Predicting pathological tumor size in prostate cancer based on multiparametric prostate magnetic resonance imaging and preoperative findings. *J Urol* 2021;205:444-451.
46. Boutros PC, Fraser M, Harding NJ, *et al.* Spatial genomic heterogeneity within localized, multifocal prostate cancer. *Nat Genet* 2015;47:736-745.
47. Salami SS, Kaplan JB, Nallandhighal S, *et al.* Biologic significance of magnetic resonance imaging invisibility in localized prostate cancer. *JCO Precision Oncol* 2019; DOI: 10.1200/PO.19.00054.
48. Valerio M, Ahmed HU, Emberton M, *et al.* The role of focal therapy in the management of localised prostate cancer: A systematic review. *Eur Urol* 2014;66:732-751.
49. Truong M, Hollenberg G, Weinberg E, *et al.* Impact of Gleason subtype on prostate cancer detection using multiparametric magnetic resonance imaging: correlation with final histopathology. *J Urol* 2017;198:316-21.
50. Li P, You S, Nguyen C, *et al.* Genes involved in prostate cancer progression determine MRI visibility. *Theranostics* 2018;8(7):1752-1765.
51. Otti, VC, Miller C, Powell RJ, *et al.* The diagnostic accuracy of multiparametric magnetic resonance imaging before biopsy in the detection of prostate cancer. *BJU Int.* 2019;123:82-90.
52. Rosenkrantz AB, Ginocchio LA, Cornfeld D, *et al.* Interobserver reproducibility of the PI-RADS version 2 lexicon: A multicenter study of six experienced prostate radiologists. *Radiology* 2016;280(3):793-804.
53. Oishi M, Shin T, Ohe C *et al.* Which patients with negative magnetic resonance imaging can safely avoid biopsy for prostate cancer? *J Urol* 2019;201:268-277.
54. Wang RS, Kim EH, Vetter JM *et al.* Determination of the role of negative magnetic resonance imaging of the prostate in clinical practice: is biopsy still necessary? *Urology* 2017;102:190-197.
55. Wysock JS, Mendhiratta N, Zattoni F, *et al.* Predictive value of negative 3T multiparametric magnetic resonance imaging of the prostate on 12-core biopsy results. *BJU Int* 2016;118:515-520.
56. Ahmed HU, El-Shater Bosaily A, Brown LC, *et al.* Diagnostic accuracy of multiparametric MRI and TRUS biopsy in prostate cancer (PROMIS): a paired validating confirmatory study. *Lancet* 2017;389:815-822.
57. Borofsky S, George AK, Gaur S, *et al.* What are we missing? False-negative cancers at multiparametric MR imaging of the prostate. *Radiology* 2018;286(1):186-195.
58. Elkhoury F, Felker E, Kwan L, *et al.* Comparison of targeted vs systematic prostate biopsy in men who are biopsy naïve: the Prospective Assessment of Image Registration in the Diagnosis of Prostate Cancer (PAIREDCAP) Study. *JAMA Surg* 2019;154(9):811-818.
59. Hogan D, Yao HHL, Kanagarajah A, *et al.* Can multi-parametric resonance imaging and prostate-specific antigen density accurately stratify patients prior to prostate biopsy? *J Clin Urol* 2022; doi.org/10.1177/20514158221084820.
60. Simmons LAM, Kanthabalan A, Arya M, *et al.* The PICTURE study: diagnostic accuracy of multiparametric MRI in men requiring a repeat prostate biopsy. *Br J Cancer* 2017;116: 1159-65.
61. Van Riel L, Jager A, Meijer D, *et al.* Predictors of clinically significant prostate cancer in biopsy-naïve and prior negative biopsy men with a negative prostate MRI: improving MRI-based screening with a novel risk calculator. *Ther Adv Urol* 2022;14:1-10.
62. Mortezaei A, Märzendorfer O, Donati OF, *et al.* Diagnostic accuracy of multiparametric magnetic resonance imaging and fusion guided targeted biopsy evaluated by transperineal template saturation prostate biopsy for the detection and characterization of prostate cancer. *J Urol* 2018;200:309-318.
63. Chung DY, Koh DH, Goh HJ, *et al.* Clinical significance and predictors of oncologic outcomes after radical prostatectomy for invisible prostate cancer on multiparametric MRI. *BMC Cancer* 2018; 18: 1057.
64. Kinnaird A, Sharma V, Chuang R, *et al.* Risk of prostate cancer after a negative MRI-guided biopsy. *J Urol* 2020; DOI: 10.1097/JU.0000000000001232
65. Moldovan PC, Van den Broeck T, Sylvester R, *et al.* What is the negative predictive value of multiparametric magnetic resonance imaging in excluding prostate cancer at biopsy? A systematic review and meta-analysis from the European Association of Urology Prostate Cancer Guidelines Panel. *Eur Urol* 2017;72:250-66.
66. Weaver JK, Kim EH, Vetter JM, *et al.* Presence of magnetic resonance imaging suspicious lesion predicts Gleason 7 or greater prostate cancer in biopsy-naïve patients. *Urol* 2016;88:119-124.
67. Ahdoon M, Wilbur AR, Reese SE, *et al.* MRI-targeted, systematic, and combined biopsy for prostate cancer diagnosis. *N Engl J Med* 2020;382: 917-28.
68. Grey A, Scott R, Shah B, *et al.* Multiparametric ultrasound versus multiparametric MRI to diagnose prostate cancer (CADMUS): a prospective, multi-centre, paired-cohort, confirmatory study. *Lancet Oncol* 2022;23:428-438.
69. Leung AK, Patil D, Howard DH and Filson CP. Payments and patient cost sharing for prostate biopsies according to image guidance, practice site and use of anesthesia. *Urol Practice* 2020;7:138-144.
70. Melo SA, Luecke LB, Kahlert C, *et al.* Glypican-1 identifies cancer exosomes and detects early pancreatic cancer. *Nature* 2015;523:177-182.
71. Steinbichler TB, Dudas J, Riecheltmann H and Skvortsova II. The role of exosomes in cancer metastasis. *Sem Cancer Biol* 2017;44:170-181.
72. Skog J, Wurdinger T, van Rijn S, Meijer DH, Gainche L, Sena-Esteves M *et al.* Glioblastoma microvesicles transport RNA and proteins that promote tumour growth and provide diagnostic biomarkers. *Nat Cell Biol* 2008;10:1470-1476.
73. Van der Vos KE, Balaj L, Skog J, Breakefield XO. Brain tumor microvesicles: insights into intercellular communication in the nervous system. *Cell Mol Neurobiol* 2011;31:949-959.
74. McKiernan J, Donovan MJ, Margolis E, *et al.* A prospective adaptive utility trial to validate performance of a novel urine exosome gene expression assay to predict high-grade prostate cancer in patients with prostate specific antigen 2-10 ng/mL at initial biopsy. *Eur Urol* 2018; 74: 731-738.
75. McKiernan J, Donovan MJ, O'Neill V, *et al.* A novel urine exosome gene expression assay to predict high-grade prostate cancer at initial biopsy. *JAMA Oncol* 2016;2(7):1-8.
76. Tutrone R, Donovan MJ, Torkler P, *et al.* Clinical utility of the exosome based ExoDx Prostate (IntelliScore) EPI test in men presenting for initial biopsy with a PSA 2-10 ng/mL. *PCPD* 2020;23(4):607-614.
77. Kretschmer A, Tutrone R, Alter J, *et al.* Pre-diagnosis urine exosomal RNA (ExoDx EPI score) is associated with post-prostatectomy pathology outcome. *World J Urol* 2022;doi.org/10.1007/s00345-022-03937-0.
78. Amparore D, Campi R, Checucci E, *et al.* Forecasting the future of urology practice: a comprehensive review of the recommendations by international and European associations on priority procedures during the covid-19 pandemic. *Eur Urol Focus* 2020;6:1032-1048.
79. de la Calle C, Fasulo V, Cowan JE, *et al.* Clinical Utility of biomarkers 4K Score, SelectMDX and ExoDx with MRI for the detection of high-grade prostate cancer. Presented at American Urological Association 2020. *J Urol Supplement* (4S) 2020;200:PD53-02.
80. Gorin M, Walsh P. Magnetic resonance imaging prior to first prostate biopsy—are we there yet? *Eur Urol* 2018;74:409-10.
81. Vickers A, Carlsson SV, and Cooperberg M. Routine use of magnetic resonance imaging for early detection of prostate cancer is not justified by the clinical trial evidence. *Eur Urol* 2020;78:304-306.
82. Vickers A, *et al.* How should molecular markers and magnetic resonance imaging be used in the early detection of prostate cancer? *Eur Urol* 2022;5:135-137.
83. Moul JW, Sant GR. How I Use It: The Exosome Diagnostics (EPI) prostate cancer biomarker utility in urology and primary care. *Can J Urol.* 2022 Aug;29(4):11224-11230. PMID: 35969726.

Where Science Intersects Innovation™

Bio-Techne® | R&D Systems™ Novus Biologicals™ Tocris Bioscience™ ProteinSimple™ ACD™ ExosomeDx™ Asuragen®



Contact Us

Exosome Diagnostics, Inc.
266 Second Ave., Suite 200, Waltham, MA 02451
1-(844) 396-7663
exosomedx.info@bio-technne.com

For research use or manufacturing purposes only. Trademarks and registered trademarks are the property of their respective owners.

MKT 211 v2.0
STRY0296642_EXO_BR_MRI_KMD

 **exosomed_x**
a **biotechne®** brand