

# How to order the ExoDx™ Prostate (IntelliScore) EPI-CE Test



## For physicians

You can learn about local laboratories offering the ExoDx Prostate EPI-CE Test or request an EPI-CE sample collection kit by contacting Exosome Diagnostics with details of your name, clinic, contact information and mailing address. The collection kit contains instructions on sample donation and return shipping. The test result will be sent to the physician/clinic when the result is available.

For more information, contact us at [exosomedx.info@bio-techne.com](mailto:exosomedx.info@bio-techne.com).



Once you're enrolled, you can order sample collection kits. Email us at [epi-emea@bio-techne.com](mailto:epi-emea@bio-techne.com) and include:

- Clinic/physician name
- Shipping address
- Contact information



Upon receipt of your collection kits, place the cold packs in the freezer. For sample collection, patients should use the included sample collection device and follow the instructions in the kit.



Complete and sign the Test Requisition Form enclosed in the kit with all information.

**That's it!** Please contact the courier listed in the kit for sample pickup and delivery to the lab. The test report will be sent to the ordering physician within 2-3 days after receipt of the sample.



## For laboratories

The ExoDx Prostate (IntelliScore) EPI-CE is intended to be used by professional users in third-party clinical laboratories competent in standard molecular biology techniques and methods. If you are with a laboratory currently providing services to primary care physicians and urologists, we invite you to contact us for more information about the ExoDx Prostate (IntelliScore) CE-IVD kit and how it can add value to your customers.



For more information about performing the EPI-CE in your own clinical laboratory, please contact [exosomedx.info@bio-techne.com](mailto:exosomedx.info@bio-techne.com) or visit [exosomedx.com/europe/exodx-prostate-epi-ce-ivd](http://exosomedx.com/europe/exodx-prostate-epi-ce-ivd).

## Questions?

If you're a **physician or a laboratory** and you need information regarding the EPI-CE Test or requesting a sample collection kit, please call our customer service team at

**+49 89 416 172 70**

**We're here to help!**


# Introducing the first and only exosomal molecular test that relies on genomic information to provide risk assessment for HGPCa.

The EPI-CE Test was developed to assist physicians to reduce unnecessary biopsies,\* and can be used as a risk assessment tool prior to MRI to determine if biopsy decision is warranted.

The EPI-CE test is based on the well-established ExoDx™ Prostate Test that has successfully operated from ExoDx's CLIA facility in the U.S. since 2016, and is reimbursed by several leading insurance providers. Both tests share the same specimen requirements, exoRNA isolation methodology, RT-qPCR detection of biomarkers and algorithm for risk score calculation. The EPI-CE test has been further developed for distribution and use in third-party

clinical laboratories in the EU and comprises a ready-to-use CE-IVD assay reagent kit and cloud-based CE-IVD software to facilitate risk score calculation.

The EPI-CE test was developed as a rule-out test (89% negative predictive value (NPV) and 92% sensitivity in the initial biopsy cohort) compared to the U.S. ExoDx Prostate test with a negative predictive value (NPV) of 91.3% and 92% sensitivity.



The 2019 NCCN guidelines include the ExoDx Prostate test (EPI) for early detection in men for both initial and prior negative biopsy.\*

\*Data on file, Exosome Diagnostics, 2021.